


Quick Start Guide



TRUFORMA Instrument Components



Assay Cartridge Components



Cartridge Preparation and Sample Volume Required

1. Confirm cartridge has not expired by reviewing the expiration date on the pouch label.
2. While wearing gloves, remove cartridge from pouch and place on a flat surface.
3. Using supplied pipette, dispense 150 μ L of the sample into sample port.
4. Ensure all of the sample has been absorbed into the sample port with no pooling. Do NOT overfill sample port.

Running a Test

1. Once the instrument is powered up, select "Run New Test" on the screen.
2. Enter the required patient information associated with the sample and select "Done".
3. A screen prompt will instruct user to insert a cartridge.
4. Place the sample-filled cartridge into the instrument by holding the sample port end with the cartridge label facing up and following the arrow-indicated insert direction displayed on the cartridge label.
5. Select sample type and "Run (Name of Assay) Test".
6. Confirm test and sample type are correct, then select "Run Test".
7. After the test is completed, the cartridge will eject. Dispose cartridge as biological waste.

Viewing Test Results

After the analysis is complete, the quantitative result — including an indicator if the result is above, below or within normal levels — will be displayed on the screen. Refer to the test cartridge package inserts or **MyZomedica** for reference range details.

Viewing Prior Results


Prior results can be accessed by selecting the Results button on the left menu screen. The last 100 results are viewable on the TRUFORMA instrument. To view all historical test results, visit **MyZomedica**.

Cleaning and Maintenance

1. Clean the Outer Case (Perform Weekly)
 - Wipe the outer case using a soft cloth, dampened with a mild, non-abrasive detergent or cleaning solution such as a 10 % bleach solution or 30% isopropyl alcohol solution. Do not spray any detergents or solutions directly onto the outer case.
2. Clean the Touch Screen (Perform Weekly)
 - Use a soft, lint-free cloth dampened with glass-cleaning fluid. Do not spray any solution directly onto the screen.
3. Replace the Air Filter (Perform Monthly)
 - Refer to section 8.3 of the TRUFORMA Instructions for Use to learn how to replace air filter.

MyZomedica – Viewing Results Remotely

1. Log in to **MyZomedica** account at myzomedica.com.
2. Choose the clinic from which you would like to see a result by selecting preferred clinic name. Confirm this is the clinic you want to access.
3. You can switch between clinics by selecting the user icon and/or user name in the upper right corner to open the user menu, then selecting the other clinic name under Recent.
4. Most recent test results may be viewed on the home screen. User can search for historical results by entering a keyword search in the search bar.
5. Select the test to view result details.



Get Support and Answers

myzomedica.com help ticket | customersuccess@zomedica.com | (734) 369-2555 Option 1
102021 Rev. 2